

**DECTRIS®**

*detecting the future*

---

MYTHEN

PILATUS

EIGER

# ***Les détecteurs HPC pour la crystallo sous conditions extrêmes***

***Andreas Förster – Application Scientist Crystallography  
ANF Réciproc – 2018-06-06***

DECTRIS Ltd.  
5405 Baden-Dättwil  
Switzerland  
[www.dectris.com](http://www.dectris.com)

# ***Outline***

***Technology***

***Direct detection***

***Photon counting***



***X-ray powder diffraction***

***The power of the threshold***




***The future of high-pressure crystallography***

# ***The meaning of HPC counting***

## ***Hybrid Pixels + Single Photon Counting***




# *Hybrid pixel = sensor + readout*





# Direct detection in hybrid pixels




## Sensor pixel

- Absorption of X-ray photons  
-> one  $e^-$ /hole pair per 3.6 eV
- Charge is captured by electric field

## Pixel electronics

- Counting of charge pulses

# Modular assembly of HPC detectors




Module active area:  $8 \times 4 \text{ cm}^2$   
100k pixels on PILATUS ( $172 \mu\text{m}$ )  
500k pixels on EIGER ( $75 \mu\text{m}$ )




# ***Superiority of direct detection***

## ***Direct detection***




## ***Indirect detection***




# Superiority of direct detection

## Direct detection



- Charge captured in electric field
- => **All photons captured**
- => **Signal remains in pixel**

## Indirect detection




- Secondary radiation scattered in scintillator
- => **Signal spread across pixels**
- => **Light partially lost**




# Superiority of direct detection

## Direct detection




- Charge captured in electric field
- => **All photons captured**
- => **Signal remains in pixel**




# Silicon and CdTe sensors

*High quantum efficiency across a wide energy range*




*Silicon for low energies*




*CdTe for high energies*

# Spatial resolution PILATUS3

*All signal contained with pixel*




20  $\mu\text{m}$  pencil beam, 10  $\mu\text{m}$  rastering




# *Spatial resolution PILATUS3 CdTe*

For CdTe as for silicon




10 µm pinhole at Ag K $\alpha$  wavelength






# 70S ribosome on EIGER X 16M



$a = 210 \text{ \AA}$ ,  $b = 450 \text{ \AA}$ ,  $c = 620 \text{ \AA}$   
Diffraction to  $2.3 \text{ \AA}$ , Y. Polikanov, UIC


# ***Photon counting – one by one***

*– Photons absorbed in sensor pixel*






# ***Photon counting – one by one***





- *Photons absorbed in sensor pixel*
- *Charge pulse proportional to energy*

# Photon counting – one by one



- Photons absorbed in sensor pixel
- Charge pulse proportional to energy
- Threshold to discard noise

# Photon counting – one by one




- Photons absorbed in sensor pixel
- Charge pulse proportional to energy
- Threshold to discard noise
- Signals above threshold are counted
  - **Suppression of dark signal**
  - Suppression of fluorescence
  - Suppression of electronic noise
- On-the-fly digitization in digital counter
  - **No readout noise**
  - Fast readout
  - High dynamic range

# ***Threshold at half photon energy***

***All signal counted within on pixel***

***Electronic noise + fluorescence filtered***




1 photon  
100% energy  
1 count

-

# Threshold at half photon energy


*All signal counted within on pixel*

*Electronic noise + fluorescence filtered*



1 photon  
100% energy  
1 count

-



1 photon  
60% energy  
1 count

40% energy

-



# Noise-free detection


100 ms

**No readout noise**  
**No dark signal**

*total number of counts: 0*

1 hour

*background radiation visible*





100K detector

# Negligible background for max. $I/\sigma(I)$

## Crystallography in vacuum

Germanate @ DLS I23.




# Dynamic range for strong reflections

## HPC detectors


*PILATUS @ 20 bit (1 million counts)*

*EIGER @ 32 bits (4 billion counts)*



## Scintillator-based detectors


*CCD @ 16 bits (65 thousand ADU)*



# Excellent spatial resolution


## HPC detectors

Sharp point spread function



## Scintillator-based detectors

Conversion dilutes signal



# ***MYTHEN2 – 1D detectors for PXRD***



*MYTHEN2 R  
for laboratories*



*MYTHEN2 X  
for synchrotrons*

---


*MYTHEN2 R/X 1D – 640 strips  
MYTHEN2 R/X 1K – 1280 strips*

*50  $\mu$ m strips  
up to 4 detectors per DCS*

# ***PXRD at elevated temperatures***

## ***Response of adenine to heat***

- Simultaneous coverage of large angular range*
- Collected data free of temporal inconsistencies*
- Low noise for accurate data*




*Courtesy Ivan Halasz (Zagreb)*



# PXRD at elevated temperatures

## *Response of adenine to heat*

- Mo source
- STOE INSITU HT2 oven
- STOE Stadi P (Debye-Scherrer)
- MYTHEN2 R 2K




Courtesy Ivan Halasz (Zagreb)

# PXRD at elevated temperatures

## *Response of adenine to heat*

- Mo source
- STOE INSITU HT2 oven
- STOE Stadi P (Debye-Scherrer)
- MYTHEN2 R 2K




Courtesy Ivan Halasz (Zagreb)

# PXRD at elevated temperatures

## Response of adenine to heat

- Mo source
- STOE INSITU HT2 oven
- STOE Stadi P (Debye-Scherrer)
- MYTHEN2 R 2K
- Two exposures of 12 min each




Courtesy Ivan Halasz (Zagreb)

# PXRD at elevated temperatures

## Multiple MYTHEN detectors

- Phase I to phase II upon heating
- Phase III above 380°C
- Phase IV briefly before phase III





# ***PILATUS & EIGER with Si sensors***

*PILATUS3 R  
for laboratory*










*EIGER R for  
laboratory*



*EIGER2 R 500K  
for laboratory XRD*

*Faster (and bigger) instruments for synchrotrons*


# Electrochemical oxidation of $\text{SrFeO}_{2.5}$



*Courtesy Werner Paulus (Montpellier)*




# ***Electrochemical cell at BM01***



*Swiss-Norwegian beamline (BM01) at ESRF*


*also on STOE Stadivari Mo  
with PILATUS 200K*

# Evolution of reflections in $0kl$ plane



$t = 5 \text{ days}$

# Threshold suppresses fluorescence




$$E_{\text{ex,Fe}} = 7.1 \text{ keV}$$

$$E_{\text{ex,Sr}} = 16.1 \text{ keV}$$

# Fluorescence suppression in the lab

## Small molecule with Europium

- $E_{em,Eu} = 5.8 \text{ keV}$
- Structure determination impossible with integrating detectors
- Excellent data with Cu radiation and threshold at 6.2 keV on EIGER R 4M




Courtesy Curtis Moore (UCSD)



# Fluorescence suppression in the lab

## Small molecule with Europium

- $E_{em,Eu} = 5.8$  keV
- Structure determination impossible with integrating detectors
- Excellent data with Cu radiation and threshold at 6.2 keV on EIGER R 4M




Courtesy Curtis Moore (UCSD)

# Crystallography at USCD

**Curtis Moore and Milan Gembicky**

– EIGER R 1M





# PILATUS3 CdTe detectors

*PILATUS3 R CdTe  
for laboratories*



*PILATUS3 X CdTe  
for synchrotrons*




$25 \times 28 \text{ cm}^2$

# Fluorescence suppression in the lab

## Charge density studies


- Dibenzyl diselenide
- Ag source
- $E_{em,Se} = 11.2-12.5$  keV
- Data collection with Ag source (22.2 keV)
- No fluorescence noise with threshold at 13.5 keV




# High-quality charge density

## Dibenzyl diselenide

- Isolevels at  $\pm 0.11 \text{ e } \text{\AA}^{-3}$



CCD




PILATUS3 R CdTe

*Courtesy Dietmar Stahlke (Göttingen)*


# High-quality charge density

## Paracyclophane

- Isolevels at  $\pm 0.047 \text{ e } \text{\AA}^{-3}$



CCD



PILATUS3 R CdTe


*Courtesy Dietmar Stahlke (Göttingen)*



# Crystallography at high pressure

## **1,2,3,4-tetrabromobenzene**

- *Material with mechanical response to heating*
- *Data collected in diamond anvil cell at 0.4 GPa*




*Courtesy Elena Boldyreva (Novosibirsk)*



# Better kit for better data


## Ten years of innovation

- Ag radiation for higher reciprocal space coverage
- HPC detector for higher-accuracy intensities
- CdTe sensor for high quantum efficiency



Oxford Diffraction Gemini  
Mo radiation  
Ruby CCD

*Acta Cryst.* (2018). E74, 613–619



Rigaku OD XtaLAB Synergy-S  
Ag radiation  
PILATUS3 R CdTe 300K

# Better kit for better data

## Ten years of innovation

- Ag radiation for higher reciprocal space coverage
- HPC detector for higher-accuracy intensities
- CdTe sensor for high quantum efficiency

|                                                                            | (e) Mo, ABSORB-7                                                       | (b) Ag, ABSORB-7                                                       |
|----------------------------------------------------------------------------|------------------------------------------------------------------------|------------------------------------------------------------------------|
| Crystal data                                                               |                                                                        |                                                                        |
| $a, b, c$ (Å)                                                              | 3.9431 (5), 10.7566 (18), 9.964 (2)                                    | 3.9390 (9), 10.781 (4), 9.944 (4)                                      |
| $\beta$ (°)                                                                | 100.557 (15)                                                           | 100.49 (3)                                                             |
| $V$ (Å <sup>3</sup> )                                                      | 415.47 (13)                                                            | 415.2 (2)                                                              |
| Radiation type                                                             | Mo $K\alpha$                                                           | Ag $K\alpha$ , $\lambda = 0.56087$ Å                                   |
| No. of reflections for cell measurement                                    | 514                                                                    | 748                                                                    |
| $\theta$ range (°) for cell measurement                                    | 2.8–22.4                                                               | 2.2–22.9                                                               |
| $\mu$ (mm <sup>-1</sup> )                                                  | 19.29                                                                  | 10.33                                                                  |
| Absorption correction                                                      | Gaussian                                                               | Gaussian                                                               |
| $T_{\min}, T_{\max}$                                                       | 0.361, 0.434                                                           | 0.486, 0.562                                                           |
| No. of measured, independent and observed [ $I > 2\sigma(I)$ ] reflections | 2116, 531, 313                                                         | 2445, 870, 496                                                         |
| $R_{\text{int}}$                                                           | 0.103                                                                  | 0.048                                                                  |
| $(\sin \theta/\lambda)_{\max}$ (Å <sup>-1</sup> )                          | 0.663                                                                  | 0.801                                                                  |
| Range of $h, k, l$                                                         | $h = -5 \rightarrow 5, k = -12 \rightarrow 11, l = -10 \rightarrow 10$ | $h = -5 \rightarrow 6, k = -14 \rightarrow 14, l = -11 \rightarrow 12$ |
| Refinement                                                                 |                                                                        |                                                                        |
| $R[F^2 > 2\sigma(F^2)], wR(F^2), S$                                        | 0.071, 0.169, 1.05                                                     | 0.037, 0.073, 0.93                                                     |
| No. of reflections                                                         | 531                                                                    | 870                                                                    |

Experiment time:  
 Ruby CCD: 32 hrs  
 PILATUS3: 6 hrs

# ***Summary***

***Direct detection for sharp spots***

***Photon counting for accuracy and speed***

***Adjustable threshold to cut fluorescence***

***CdTe sensors for high X-ray energies***

**DECTRIS®**

*detecting the future*

---

***Merci pour votre attention!***

***[www.dectris.com](http://www.dectris.com)***

***[andreas.foerster@dectris.com](mailto:andreas.foerster@dectris.com)***

**DECTRIS Ltd.**  
5405 Baden-Dättwil  
Switzerland  
**[www.dectris.com](http://www.dectris.com)**